

onSIGHT™ A.V.S

AMBULATORY VIDEO EEG SOLUTION



Monitor your patient's home ambulatory studies with **CONFIDENCE!**

EEG with Video recording in "Real Time"!

- > True web application for remote monitoring
- > Access through any commercial browser
- > View all data on computer or tablet
- > Office, home, or mobile access

By applying the very latest in technological innovations, Compumedics is now able to provide truly informative real-time interaction with ambulatory recordings. OnSight recordings are collected and reviewed using Compumedics Profusion EEG software, the same software used in our sophisticated laboratory systems around the world.



High Definition Camera

- > Integrated infrared light source
- > 150 degree image framing
- > 30 frames per second

Easy to transport

- > Reinforced wheeled case with retractable handle
- > Less than 25.4 lbs (11.52kg) total weight

Time Savings

- > Simple set-up: 5 minutes to enter patient data and program internal memory card
- > Less than 20 minutes download time (USB3) for typical 72 hour EEG/video recording (60GB)

Simple Two-Step Recording Start

1. Patient plugs the ONsight system into a standard home power outlet
2. Extend the camera mount and point the camera towards patient



Protective camera cover included



Recording Duration Options

- Full 24 hours - 4 internal lithium AA batteries
- Full 72 hours - with external battery pack
- Your patient never has to interrupt the recording to change the batteries.

Data Security

- > Two copies of recording are preserved - Internal memory card and PC
- > HIPAA 256 AES encryption
- > Viewed over secured cellular network



Patient wearing ONsight recorder with external battery pack.

profusion EEG[™]
Unlimited Potential in Routine and LTM EEG

Best in Class EEG Application Software

Compumedics Profusion EEG is the true next generation clinical and LTM software package that offers unprecedented ease of use both in the EEG lab, and in post-recording analysis and review. Developed in close consultation with EEG clinicians and neurologists, it is designed for easy navigation and editing of up to 45 days of continuously recorded EEG.

profusion nexUS 360[™]

A Revolution in Laboratory Management

Fully web-based management tool that provides seamless, automatic data archiving, study status tracking, complete program metrics and anytime, anywhere data access via standard browser interface. Profusion EEG review software is integrated into the NeXus 360 application. The perfect way to optimize your new ONsight program!

All specifications are subject to change without notice.
Check with your local Compumedics representative for product availability.
ONsight is a trademark of Compumedics. Compumedics and the Compumedics logo are all registered trademarks of Compumedics Limited.